

# Topigenix



**Physical Description:** Topigenix Wound Cream is an acidic zinc-saline cream for external wound management. It aids the body during the healing process by providing a moist, acidic wound environment. Topigenix Wound Cream provides the optimal physical environment, (moisture, pH-level) for the external management of wounds.

**Indications:** Topigenix Wound Cream is intended to manage dermal lesions such as lacerations, post-operative (surgical) wounds, partial and full-thickness wounds, 1st and 2nd degree burns, ulcers (diabetic, venous stasis, and pressure). Topigenix Cream can be used in conjunction with dressings that absorb fluid (i.e. gauze, gel, alginate, foam and hydrocolloid).

**Directions:** Apply to affected area of the skin as recommended below or as directed by the physician.

**For Acute Wounds:** Apply to the tissues surrounding the wound. With wounds such as incisions or puncture wounds, it is recommended that the wound first be cleaned with Topigenix Skin & Wound Cleanser and that the wound cavity or incision be filled with Topigenix Wound Gel. Topigenix Cream should be applied to the wound margins to provide nutrients promoting healing. The wound should then be covered with a dressing. Once the wound has closed and healing has begun, the wound cream should be applied directly over the wound with continued application for up to 28 days to promote continued healing and tissue regeneration.

## Topigenix Wound Cream

1 oz Tube (28.35 gm)

- Provides a moist and soothing healing environment
- Provides essential elements for healing
- Provides a low pH
- Non-toxic, non-sensitizing, non-irritating

**For Chronic Wounds:** When treating chronic wounds, such as pressure ulcers, the open wounds and surrounding tissues should be cleansed with Topigenix Skin & Wound Cleanser and the wound bed filled with Topigenix Wound Gel. Topigenix Cream should be applied to the surrounding wound margins to provide additional nutrients to these tissues and to prevent wound expansion. The wound should then be covered with a dressing. This treatment process should be repeated every other day or as directed. After the wound bed has closed (usually within 5-7 days of treatment) the Wound Cream should be applied 2-3 mm thick directly over the wound once daily. Continued application to the wound for up to 28 days after therapy will help to promote tissue remodeling and reduce scarring.

**Contraindications:** Topigenix Wound Cream is contraindicated in the following cases:

- Known or suspected allergic reactions to any of the ingredients
- Third degree burns
- Any wound where the dermis is severely damaged or is missing.

**Warning:** For external use only. Do not use in the eyes.

**Shelf Life and Storage:** Shelf Life of 36 months (3 years). Use within 3 months of opening. Store at controlled room temperature (41-86° F).

Manufactured in the USA exclusively for:  
Healthcare Industry Group, Inc., Naples, FL  
[www.healthcareindustrygroup.com](http://www.healthcareindustrygroup.com)  
[Topigenix@healthcareindustrygroup.com](mailto:Topigenix@healthcareindustrygroup.com)  
Phone: 630-450-0639

**HIG** Healthcare Industry Group